

7000



Fully powered, hydraulic
examination chair/
minor **procedures table**

Reliance®

7000



Flexibility for examinations or minor procedures.

In addition to its optional pillow headrest, the Model 7000 reclines to the near flat position, which enables it to perform double duty as an examination chair and a minor procedures table. A surgical headrest is optional.

Powered Base and Articulating Top

Use finger switches and corded foot switch to raise, lower or recline chair into any position without breaking contact with patient. Eliminates wasted movements. An auto return button automatically returns the chair to the low, upright position.

Memory Positioning/Safety Switch

A memory switch allows you to instantly pre-select the chair's most used position. The patented programmable safety switch disables side and foot switches which prevents children from operating the chair.

Membrane Switches

Membrane switches enhance the Model 7000's sleek design and contribute to its durability and safety. They last 10-times longer than mechanical switches, eliminate the need for high voltage and allow limited programming of chair functions. Back lighting eases work in a dimmed room. Mechanical switches available as option.

Dual Rotation Locks

Easily accessible levers located on both sides of seat enable practitioner to swivel chair 330° without changing position.

Large Hinged Footstep

The 8" x 13" footstep supports patient's legs, and provides safe ingress and egress. Tilts out of way. Padded underneath for comfort.

Headrest Choices

Standard concave upholstered headrest offers superior support and a wide range of positions. Locks firmly with a single hand. Options include a non-articulating pillow headrest and an articulating surgical headrest.

Armrests

Armrests lock in place. Release mechanism allows convenient ingress and egress from either side of chair. Padded underneath for comfort.

Upholstery

Match office or treatment room décor with our selection of 14 easy-to-maintain vinyl colors. Other colors can be supplied at additional cost. Available with either smooth or the designer-style cushions shown on cover. Base is available in #20 Pearl. All materials are fire retardant. To request a color selection card, contact your dealer or call Reliance at 1-800-735-0357.

The CSA Mark of Quality

The Model 7000 Chair has passed rigorous testing by CSA Laboratories, which assures product quality and safety for you and your patients. The Model 7000 chair carries the CSA, C/US Certificate of Compliance.

Dimensions and Specifications

Base dimensions: 22 3/4" x 24 3/4". Overall width: 24 1/4". Overall depth in chair position (excluding headrest): 36". Table position: 60". Seat: 18 1/2" (wide) x 17" (deep). Distance between armrests: 19". Choice of two base heights: HIGH: 22 3/4" lowered, 34" raised; LOW: 20" lowered, 29" raised. Shipping weight: 325 lbs. Electrical: (Standard) 115 volts, 60 Hz, 8.0 amperes total; (Optional) 230 volts, 50 Hz, 4.0 amperes.

