

150



Affordable,
small footprint
ENT Cabinet

Reliance®



Spacious Work Surface

Laminate work surface cleans up quickly and provides the practitioner with ample room for charting and maneuvering.

Convenient Work Center

Instruments are strategically placed above the work surface but within practitioner's reach. This provides a large work surface uncluttered by bottles and instruments.

Suction and Pressure

1000 cc disposable suction container is equipped with an automatic shut-off to prevent spills. Pressure hose comes with a thumb piece regulator. Holders on side of cabinet store hoses out of way.

Control Console

On/off circuit breaker switch, vacuum and pressure pump switches, and LED instrument charging indicators are available on either the right or left side for the practitioner's convenience.

Drawers

Drawer inserts keep instruments segregated and organized.

Gast Pumps

Heavy duty Gast pressure and suction pumps feature a five micron air filter and air pressure regulator that eliminates particulates. Optional foot pedal controls pumps.

Instrumentation and Options

Standard instrumentation includes two 3.5 volt rechargeable instruments, including your choice of heads, and three DeVilbiss glassware sets with two stainless steel cotton containers. Options include corded footswitch, lift-up shelf and corded instruments.

Dimensions and Specifications

Overall height: 48 3/8". Height to work surface: 33 3/8". Width: 24 3/8". Depth: 20 1/4". Shipping weight: 260 lbs.
Electrical: (Standard) 115 volts; (Optional) 230 volts.



Precision by Tradition