

400

Treatment cabinet
offers **easy access**
to instruments



Reliance®

Shown with optional corded instruments

400

Classic treatment cabinet.

The Model 400 provides easy accessibility to instruments and generous drawer space.



Spacious Work Surface

Laminate surface cleans up quickly and provides the practitioner with ample room for charting and maneuvering.

Convenient Work Center

Instruments are strategically placed within practitioner's reach. Work center also features two built-in countertop rechargeable wells.

Suction and Pressure

800 cc disposable suction container is equipped with an automatic shut-off to prevent spills. Pressure hose is equipped with thumb piece regulator. Holders on side of cabinet store hoses out of way.

Control Console

On/off circuit breaker switch, vacuum and pressure pump switches, and LED instrument charging indicators are available on either the right or left side for the practitioner's convenience.

Full-Extension Drawers

Generously proportioned removable full-extension drawer inserts are divided to keep instruments segregated.

Warming Drawer

Select the Reliance concealed mirror warmer or the Premier beaded warmer to prevent condensation on instruments.

Gast Pumps

Heavy duty Gast pressure and suction pumps feature a five micron air filter and air pressure regulator that eliminates particulates. Optional foot pedal controls pumps.

Choice of Laminate

Choice of four laminate colors may be matched with an accent color of black or white. Laminates add durability and simplify cleaning.

Instrumentation and Options

Standard instrumentation includes two Welch Allyn 3.5 volt rechargeable instruments and a complete DeVilbiss glassware set with two stainless steel cotton containers. Options include foot switch, chart holder, concealed mirror warmer, and debris tray and holder.

Dimensions and Specifications

Overall height: 44". Height to work surface: 31". Width: 29". Depth: 18". Shipping weight: 150 lbs. Electrical: (Standard) 115 volts; (Optional) 230 volts.



Precision by Tradition