

425



Treatment Cabinet

Reliance®

Shown with optional dual light source and bead warmer

425

Treatment Cabinet

The Model 425 cabinet's generously proportioned full-extension drawers provide ample storage for instruments and equipment.



Spacious Work Surface

Seamless surface cleans up quickly and provides the practitioner with ample room for charting and maneuvering.

Elevated Work Center

Elevated tower provides easy access to instruments and prevents them from cluttering the work surface and coming in contact with contaminants. Work center also provides area for dispensing gloves and tissues.

Suction and Pressure

1000 cc disposable suction container is located on a swivel arm for easy changing. Automatic shutoff prevents spills. Pressure hose is equipped with thumb piece regulator. Holders on side of cabinet store hoses out of way.

Control Console

On/off switch, vacuum and pressure pump switches, and LED instrument charging indicators are kept at the practitioner's fingertips.

Full-Extension Drawers

Generously proportioned full-extension drawers are equipped with removable plastic liners to keep instruments segregated and to facilitate cleaning.

Warming Drawer

Select the Reliance concealed mirror warmer or the Premier beaded warmer to prevent condensation on instruments.

Gast Pumps

Gast pressure and suction pumps, housed in a soundproofed compartment at the bottom of the unit, feature a five micron air filter to eliminate particulates. Optional foot pedal controls pumps.

Choice of Laminate on Drawer Fronts

Match office décor with a choice of 10 laminate colors from the Wilson Art Nebula series. Laminates add durability and simplify cleaning.

Instrumentation and Options

Standard instrumentation includes two Welch Allyn 3.5 volt rechargeable instruments and a complete DeVilbiss glassware set with two stainless steel cotton containers. Options include foot switch, Reliance concealed drawer warmer or the Premier beaded warmer, a debris tray and cover, debris tray holder, and a built-in port for fiber optic light source.

Dimensions and Specifications

Overall height: 45". Width: 32". Depth: 20". Height to work surface: 31". Drawer sizes inside: 3@ 13 1/2" x 15" x 2 3/8" and 1@ 18 1/2" x 15" x 8 1/2" Shipping weight: 350 lbs. Electrical: 110 volts, 60 Hz, 9 amperes (12 amperes with light source). Patent 5,399,077. UL Listed.



Precision by Tradition