

### Digital cameras evaluated by Haag-Streit on BP 900 (Photo adapter DP 01)

Digital camera model	Recommended by Haag-Streit	Possible	Not possible	Step down ring	Haag-Streit intermediate ring	Remote release	AC Adapter	Remarks
Canon Power Shot A640	X			Canon LA-DC58F	1008548	not possible	-	-
Canon Power Shot A630	X			Canon LA-DC58F	1008548	not possible	-	-
Canon Power Shot S80	X			Canon LA-DC20	1006974	not possible	ACK-DC20	-
Canon Power Shot S60	X			Canon LA-DC10	1006974	Canon WL-DC100	ACK-700	remote release recommended
Canon Power Shot S70	X			Canon LA-DC10	1006974	Canon WL-DC100	ACK-700	remote release recommended
Canon Power Shot A70	X			Canon LA-DC52C	1006829	not possible	ACK-600	-
Canon Power Shot A80	X			Canon LA-DC52D	1006829	not possible	ACK-600	-
Canon Power Shot A510	X			Canon LA-DC52F	1006829	not possible	ACK-800	-
Canon Power Shot A520	X			Canon LA-DC52F	1006829	not possible	ACK-800	-
Nikon Coolpix 4300	X			Nikon UR-E4	1006767	Nikon MC-EU1	EH-53	-
Nikon Coolpix 4500	X			not needed	1006767	Nikon MC-EU1	EH-53	-
Sony Cybershot DSC-W1	X			Sony VAD-WA	1007138	not possible	AC-LS5	-
Sony Cybershot DSC-W5	X			Sony VAD-WA	1007138	not possible	-	-
Nikon Coolpix 8400		X		Nikon UR-E16	1006767	ML-L3 (IR)	EH 54	may be inconvenient for head; remote release not recommended
Sony DSC-V3		X		Sony VAD-VHA	1006886	not possible	AC-LS5	slight vignetting with optic zoom
Sony Cybershot DSC-W7		X		Sony VAD-WA	1007138	not possible	-	Image quality not checked
Canon PowerShot A 700			X	Canon LA-DC 58G	1006886	not possible	-	-
Nikon Coolpix 5700			X	-	-	-	-	-
Sony DSC-V1			X	-	-	-	-	no step down ring