

### Digital cameras evaluated by Haag-Streit on BX 900

Digital camera model	Recommended by Haag-Streit	Possible	Not possible	Haag-Streit intermediate ring	Remote release	AC Adapter	Remarks
Canon EOS 30D	X			3000827	1007306	ACK-E2	
Canon EOS 10D	X			3000827	1007306	ACK-E2	replaced by 20D
Canon EOS 20D	X			3000827	1007306	ACK-E2	replaced by 30D
Canon EOS D30	X			3000827	1007306	-	USB1 only
Canon EOS D60	X			3000827	1007306	-	USB1 only
Nikon D200	X			3000828	1007305	EH-6	-
Nikon D1x	X			3000828	1007305	-	Fire wire connector inconveniently placed on the backside
Nikon D2x	X			3000828	1007305	-	-

Other camera models might be adapted to the system but this leads to increased costs and delivery time because of adaptations necessary and image quality can not be granted. Please contact your local HS distributor to determine the camera that fits your needs best.